



EMPAT™ QC SAMPLER

For Automatic Aseptic Sampling in cGMP Production

MINIMIZE CONTAMINATION RISK

- On-line aseptic QC sampling during batch processing
- Pre assembled single use set with 7-sample bags and aseptic connector

AUTOMATE CLEAN ROOM OPERATION

- Day-to-day reliable performance, following customer sampling SOP
- Walk away operation, process and sample data monitored and archived
- Critical process parameters and events alarmed, and archived

ADDITIONAL FEATURES

- Input signals from auxiliary equipment can be archived
- 21CFR compliant software includes user log-on, multi-level password protection and e-signature

ENVIRONMENTAL MONITORING

- Registration and (bar code) labeling of environmental monitoring samples collected during batch processing
- On-board storage enclosure for EM samples can be made available (optional)





Image for illustrative purposes only